THUASNE®

ROM-R

fr	Attelle postopératoire avec limitation de flexion/extension.....	4
en	Postoperative splint with flexion/extension limitation.....	4
de	Postoperative Knieschiene mit Beugeoder Streckbegrenzung.....	4
nl	Postoperatieve brace met buig-/ strekbepierking.....	5
it	Tutore post-operatorio con limitazione di flessione/estensione.....	5
es	Férula postoperatoria con limitación de la flexión / extensión.....	5
pt	Tala pós-operatória com limitação de flexão/extensão.....	6
da	Postoperatorisk skinne med begrænsning af fleksion/ekstension.....	6
fi	Leikkauksen jälkeinen tuki koukistuksen/ojennuksen rajoittamiseksi.....	6
sv	Postoperativ ortos med flexions-/extensionsbegränsning.....	7
el	Μετεγχειρητικός νάρθηκας με περιορισμό της κάμψης / έκτασης.....	7
cs	Pooperační ortéza s nastavením rozsahu flexe/extenze.....	8
pl	Orteza pozbiegowa z ograniczeniem kąta zgięcia/wyprostu.....	8
lv	Pēcooperācijas šina ar ceļa saliekšanas un atliekšanas ierobežošanu.....	8
lt	Lenkimą / ištiesimą ribojantis pooperacinis įtvaras.....	9
et	Painutust-sirutust piirav operatsioonijärgne lahas.....	9
sl	Pooperacijska opornica z omejitvijo fleksije/ekstenzije.....	9
sk	Pooperačná dlaha s obmedzením flexie/extenzie.....	10
hu	Műtét utáni sín hajlítás/nyújtás korlátozással.....	10
bg	Следоперативна шина с ограничение на сгъване/разгъване.....	10
ro	Dispozitiv/Atelă postoperatorie cu reglarea flexiei/extensiei.....	11
hr	Postoperativna ortoza s ograničenjem fleksije/ekstenzije.....	11
zh	术后限制弯曲/伸展的夹板.....	11
ar	جبيرة ما بعد الجراحة ذات مفصل قابل لتعديل الزاوية.....	12

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso - Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Råd vid användning - Οδηγίες χρήσης - Pokyny k používání - Szczegółowe użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare - Инструкция по эксплуатации - Upute za uporabu - 使用指南 - إرشادات الاستعمال

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






© Thuasne 2017002 (2022-12)
©Studio Caterin

THUASNE
120, rue Marius Auphan
92300 Levallois-Perret
France

www.thuasne.com

fr	Taille unique
en	Universal size
de	Einheitsgröße
nl	Eén maat
it	Taglia unica
es	Talla única
pt	Tamanho único
da	Kun én størrelse
fi	Yksi koko
sv	En storlek
el	Ένα μέγεθος
cs	Univerzální velikost
pl	Rozmiar uniwersalny
lv	Universāls izmērs
lt	Vienas dydis
et	Üks suurus
sl	Enotna velikost
sk	Univerzálna veľkosť
hu	Egy méret
bg	Един размер
ro	Mărime universală
hr	Univerzalna veličina
zh	均碼
ar	مقاس وحيد

fr	Immobilisation	Contrôle du mouvement
en	Immobilisation	Motion control
de	Ruhigstellung	Bewegungsführung
nl	Immobilisatie	Bewegingscontrole
it	Immobilizzazione	Controllo del movimento
es	Inmovilización	Control del movimiento
pt	Imobilização	Controlo do movimento
da	Støtte	Kontrol af bevægelsen
fi	Immobilisoi	Liikkeen hallinta
sv	Immobilisering	Rörelsekontroll
el	Ακινητοποίηση	Έλεγχος της κίνησης
cs	Znehybnění	Kontrola pohybu
pl	Unieruchomienie	Kontrola ruchu
lv	Imobilizācija	Kustības kontrole
lt	Imobilizavimas	Judesio kontrolė
et	Liikumatuks muutmine	Liikuvuse kontroll
sl	Imobilizacija	Nadzor gibanja
sk	Imobilizácia	Vymedzenie rozsahu pohybu
hu	Rögzítés	Mozgásszabályozás
bg	Обездвижаване	Контрол на движението
ro	Imobilizare	Controlul mișcării
hr	Imobilizacija	Kontrola pokreta
zh	固定	运动控制
ar	التثبيت	التحكم في الحركة

fr

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en

POSTOPERATIVE SPLINT WITH FLEXION/EXTENSION LIMITATION

Description

Destination : the splint is aimed at patients with a severe ligament injury or in postoperative treatment. The patient's thigh circumference must correspond to the fitting range (up to 80 cm – 31.5"), for a minimum leg length of 58 cm (22").

The splint consists of :

- 2 height-adjustable (without tools) rigid side uprights (aluminium and polypropylene) **(A)** ;
- an adjustable monocentric hinge on each upright (limitation of movements and flexion/extension blocking) **(B)** ;
- 4 support straps (polyamide) with buckles (polyoxymethylene) **(C)** ;
- condylar inserts and removable foam padding (ethylene-vinyl acetate, polyamide, polyester) for patient comfort **(D)**.

The product is totally open and easy to put on, with buckles for quick, easy fitting and removal.

The product can be worn on either leg. It is available in one size only.

Properties/Mechanism of action

- Stabilisation of knee joint ligaments provided by the side uprights.
- Adjustment of splint flexion (0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120°) or extension (0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°).
- Total immobilisation possible at 0°, 15°, 30°, 45°, 60°, 75° or 90°.

Indications

- Immobilisation of the knee during the postoperative and/or functional rehabilitation period (including ligament, patellar tendon, joint cartilage and meniscus injuries, as well as fixed or stable fractures of the tibial plateau, condyles, proximal tibia or distal femur).
- Immobilisation of the knee following an injury.
- Conservative treatment of knee ligament injuries and/or ruptures.
- Knee instability.

Contraindications

- Unstable or surgical fractures.
- Fractures outside the product's effective zone (proximal femur fracture, distal tibia or fibula fracture).
- Known allergy to any of the components.
- Major venous thromboembolic history without thromboprophylaxis.

Undesirable side-effect

- Wearing this product can cause skin reactions (redness, erythema, etc.) in some people.
- Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Precautions

- Closely follow your healthcare professional's prescription and recommendations for use.
- Verify the product's integrity before every use.
- Do not place the product in direct contact with broken skin. Carefully cover any skin lesions with a suitable dressing.
- It is recommended that the venous thromboembolic risk be assessed before any immobilisation is indicated. Follow the advice of a healthcare professional.
- In the event of discomfort, significant restriction or unusual sensations, remove the product and seek the advice of a healthcare professional.
- In the event of a change in product performance, contact the healthcare professional who prescribed or supplied the product.
- This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.
- For hygiene and performance reasons, do not re-use the product for another patient.
- Remove the product before taking a bath or shower.
- Remove the product before any imaging exams.
- Store at room temperature, preferably in its original packaging.

Instructions for use/Application

The product should be adjusted and fitted for the first time by a qualified healthcare professional.

- Position the open splint under the leg. Line up the centre of the hinge with the centre of the knee (figure 1).
- Shaping the uprights: the uprights can be shaped if necessary to prevent any pressure points or splint migration.
 - Make sure the splint length is set to the minimum (the splint is telescoped), then place it on a flat surface to protect the hinge.
 - Apply gentle pressure in the desired direction. Repeat the process on the other upright if necessary.
- Adjusting the height:
 - The splint height can be adjusted from 48 cm to 66 cm (19" – 26").
 - Unlock the end of the upright and slide to adjust to the desired length (figure 2).
 - Apply the same length for the uprights on both sides.
 - It is possible to lengthen the splint at the thigh and/or calf.
- Fastening the straps:
 - Engage the buckles (figure 3.1).
 - Tighten the straps in the correct order (figure 3.2).
 - To adjust the strap length, remove the touch fastener, cut the strap using scissors then replace the touch fastener.
 - Optional: to secure the straps more firmly over the hinges, additional touch fasteners are available.
- Adjusting the hinge:
 - Release the blue catch identified by the .
 - The leg can be immobilised from 0° to 90° by increments of 15°. To adjust the immobilisation angle, position the product at the desired inclination and engage the blue catch again (figure 4).
 - The extension angle can be adjusted from 0° to 105° by increments of 15°. The flexion angle can be adjusted from 0° to 120° by increments of 15°. To set the authorised range of motion, push the black buttons located on each side and slide them to the desired angle (figure 5).
 - The setting must be exactly the same on both hinges.

f. Removing the splint:

- To remove the splint, release the buckles: raise the buckle and apply a rotating movement to release it (figure 6).

g. Optional – Drawstrings are available to avoid any need for manipulation by the patient. Insert the drawstrings through each button, tighten and cut off the excess.

Care/Maintenance

When necessary, remove the foam padding, hand-wash using mild soap, squeeze out excess water and air dry away from direct heat sources.

Do not use detergents, fabric softeners or aggressive products (products containing chlorine).

The uprights may be cleaned with a damp cloth. Do not immerse the product when washing.

The product should be disposed of as per local regulations.

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